

Morinaga Food Allergen Lateral Flow
Rapid Test Pro II for Egg

Rev. May. 2025

Lot number from #PEG036 -

Group	Ingredient	Test Result	Group	Ingredient	Test Result	Group	Ingredient	Test Result
Cereals	Rice powder	-	Seeds	Peanut	-	Fruits	Orange	-
	Polished rice	-		Roasted almond	-		Kiwi	-
	Sprouted brown rice	-		Roasted cashew nut	-		Apple	-
	Glutinous rice	-		Roasted macadamia nut	-		Peach	-
	Buckwheat	-		Roasted pistachio nut	-		Banana	-
	Japanese barnyard millet	-		Roasted walnut	-	Coffee・Cocoa・Tea	Coconut milk	-
	Foxtail millet	-		Roasted pecan nut	-		Coffee bean	-
	Common millet	+		White sesame	-		Cocoa powder	-
	Wheat	-	Sea foods	Black sesame	-		Green tea	-
	Rye	-		Salmon	-	Seaweeds	Roasted laver	-
	Rye	-		Mackerel	-		Hijiki (boiled)	-
	Barley	-		Squid	-		Brown seaweed (raw)	- *
	Malt	-		Octopus	-		Kelp 200-fold extract	-
	Oat	-		Black tiger prawn	-	Spices	Cumin	-
	Adlay	-		Snow crab	-		Poppy seeds	-
	Corn	-		Scallop adductor	-		Black pepper	-
	Corn flour	-		Pollack roe	-		White pepper	-
	Amaranth	-		Salmon roe	-		Red pepper	-
	Quinoa	-		Beef	-		Ginger	-
Tubers・Starch	Sago palm powder	-	Meats	Chicken	-		Cinnamon	-
	Yam	-		Pork	-		Coriander	-
	Potato starch	-		Chicken egg	+		Turmeric	-
Beans	Soybean	-	Eggs	Processed chicken egg	+	Polysaccharide Thickener	Low methoxyl pectin 50-fold extract	-
	Red bean (azuki)	-		Dried chicken egg	+		Guar gum 1000-fold extract	-
	Common bean	-	Dairy products	Milk	-	Others	Gelatin	-
	Concord pole	-		Skim milk	-			
	Broad bean	-		Tomato	-			
			Green stuff	Onion	-			
				Spinach	-			
			Fungus	Mushroom	-			

Test Result	
-	Negative
+	Positive

*The Test Solution may be difficult to flow on a test stick.

These samples were tested by following the instruction of Rapid Test Pro II (using Recommended Extraction Method), except Kelp 200-fold extract,Polysaccharide and Thickener samples.

These results can be changed by the difference of kit's lot number or the sample's species or production area.