

Morinaga Food Allergen ELISA
Total Milk ELISA Kit II

Rev. July. 2024

Lot number from #014- (the last three digits)

Group	Ingredient	Value μg/g	Group	Ingredient	Value μg/g	Group	Ingredient	Value μg/g
Cereals	Rice powder	<0.31	Sea foods	Salmon	<0.31	Coffee・ Cocoa・Tea	Coffee bean	<0.31
	Polished rice	<0.31		Mackerel	<0.31		Cocoa powder	<0.31
	Sprouted brown rice	<0.31		Tuna	<0.31		Green tea	<0.31
	Red rice	<0.31		Dried bonito	<0.31		Red tea	<0.31
	Black rice	<0.31		Squid	<0.31	Seaweeds	Green laver	<0.31
	Glutinous rice	<0.31		Octopus	<0.31		Roasted laver	<0.31
	Buckwheat	<0.31		Black tiger prawn	<0.31		Brown seaweed	<0.31
	Japanese barnyard millet	<0.31		Snow crab	<0.31		Dried kelp	<0.31
	Foxtail millet	<0.31		Japanese spiny lobster	<0.31	Spices	Rosemary	<0.31
	Common millet	<0.31		Lobster	<0.31		Vanilla beans	<0.31
	Wheat	<0.31		Clam	<0.31		Cumin	<0.31
	Rye	<0.31		Abalone	<0.31		Poppy seeds	<0.31
	Barley	<0.31		Scallop adductor	<0.31		Fennel	<0.31
	Malt	<0.31		Pollack roe	<0.31		Fapaniese peper	<0.31
	Oat	<0.31		Salmon roe	<0.31		Black pepper	<0.31
	Adlay	<0.31	Meats	Raw beef	<0.31		White pepper	<0.31
	Corn	<0.31		Roasted beef	<0.31		Cardamon	<0.31
	Corn flour	<0.31		Beef liver	<0.31		Sage	<0.31
	Amaranth	<0.31		Raw chicken	<0.31		Clove	<0.31
	Quinoa	<0.31		Roasted chicken	<0.31		Red pepper	<0.31
Tubers・ Starch	Sago palm powder	<0.31		Chicken liver	<0.31		Horseradish	<0.31
	Potato	<0.31		Raw pork	<0.31		Wasabi (green horseradish)	<0.31
	Potato flake	<0.31		Roasted pork	<0.31		Ginger	<0.31
	Yam	<0.31		Pork liver	<0.31		Curry powder (spice mix)	<0.31
Beans	Soybean	<0.31	Eggs	Chicken egg	<0.31		Nutmeg	<0.31
	Red bean (azuki)	<0.31		Processed chicken egg	<0.31		Cinnamon	<0.31
	Common bean	<0.31		Dried chicken egg	<0.31		Coriander	<0.31
	Green pea	<0.31	Dairy products	Milk			Basil	<0.31
	Runner bean	<0.31		Butter	>20		Thyme	<0.31
	Kidney bean	<0.31		Skim milk			Laurel	<0.31
	Concord pole	<0.31	Green stuff	Tomato	<0.31		Dried tangerine peel	<0.31
	Black bean	<0.31		Onion	<0.31		Turmeric	<0.31
	Broad bean	<0.31		Spinach	<0.31	Polysaccharide Thickener	Carrageenan (i-type) 100-fold extract	<1.56
Seeds	Peanut	<0.31	Fungus	Garlic	<0.31		Carrageenan (k-type) 50-fold extract	<0.78
	Roasted almond	<0.31		Shiitake	<0.31		Carrageenan (λ-type) 50-fold extract	1.96
	Roasted pecan nut	<0.31		Matsutake	<0.31		Gellan gum 50-fold extract	<0.78
	Roasted cashew nut	<0.31		Mushroom	<0.31		Locust bean gum 200-fold extract	<3.12
	Roasted macadamia nut	<0.31	Fruits	Orange	<0.31		Low methoxyl pectin 50-fold extract	<0.78
	Roasted pistachio nut	<0.31		Kiwi	<0.31		Xanthane gum 50-fold extract	<0.78
	Roasted hazel nut	<0.31		Apple	<0.31		Guar gum 1000-fold extract	<15.6
	Roasted walnut	<0.31		Peach	<0.31	Others	Gelatin	<0.31
	White sesame	<0.31		Banana	<0.31		Yeast extract	<0.31
	Roasted white sesame	<0.31		Raisin	<0.31			
	Black sesame	<0.31		Coconut powder	<0.31			
	Roasted black sesame	<0.31		Coconut milk	<0.31			
	Ginkgo nut	<0.31						
	Pine nut	<0.31						
	Matrimony vine seed	<0.31						
	Gardenia seed	<0.31						
	Hemp seed	<0.31						

N.T. : Not Tested

Reaction positive

※These samples were measured by following the instruction of this ELISA Kit.

※These results can be changed by the difference of kit's lot number or the sample's species or production area.