

Morinaga Food Allergen Lateral Flow    Food reactivity data  
Rapid Test ProⅡ for Egg

Rev. May. 2025

Lot number from #PEG036 -

Group	Ingredient	Test Result	Group	Ingredient	Test Result	Group	Ingredient	Test Result	
Cereals	Rice powder	-	Seeds	Peanut	-	Sea foods	Salmon	-	
	Polished rice	-		Roasted almond	-		Mackerel	-	
	Sprouted brown rice	-		Roasted cashew nut	-		Squid	-	
	Glutinous rice	-		Roasted macadamia nut	-		Octopus	-	
	Buckwheat	-		Roasted pistachio nut	-		Black tiger prawn	-	
	Japanese barnyard millet	-		Roasted walnut	-		Snow crab	-	
	Foxtail millet	-		Roasted pecan nut	-		Scallop adductor	-	
	Common millet	+		White sesame	-		Pollack roe	-	
	Wheat	-		Black sesame	-		Salmon roe	-	
	Rye	-	Green stuff	Corn	-	Meats	Beef	-	
	Barley	-		Tomato	-		Chicken	-	
	Malt	-		Onion	-		Pork	-	
	Oat	-		Spinach	-	Eggs	Chicken egg	+	
	Adlay	-	Orange	-	Processed chicken egg		+		
	Corn flour	-	Kiwi	-	Dried chicken egg		+		
	Amaranth	-	Apple	-	Dairy products		Milk	-	
	Quinoa	-	Peach	-		Skim milk	-		
Tubers・Starch	Sago palm powder	-	Fruits	Banana	-	Coffee・Cocoa・Tea	Coffee bean	-	
	Yam	-		Coconut milk	-		Cocoa powder	-	
	Potato starch	-		Fungus	Mushroom		-	Green tea	-
Beans	Soybean	-	Seaweeds	Roasted laver	-		Spices	Cumin	-
	Red bean (azuki)	-		Hijiki (boiled)	-	Poppy seeds		-	
	Common bean	-		Brown seaweed (raw)	- *	Black pepper		-	
	Concord pole	-		Kelp 200-fold extract	-	White pepper		-	
	Broad bean	-				Red pepper		-	
								Ginger	-
								Cinnamon	-
								Coriander	-
								Turmeric	-
						Polysaccharide Thickener	Low methoxyl pectin 50-fold extract	-	
							Guar gum 1000-fold extract	-	
						Others	Gelatin	-	

Test Result	
-	Negative
+	Positive

\*The Test Solution may be difficult to flow on a test stick.

These samples were tested by following the instruction of Rapid Test ProⅡ (using Recommended Extraction Method), except for Kelp 200-fold extract,Polysaccharide and Thickener samples.

These results can be changed by the difference of kit's lot number or the sample's species or production area.