

Morinaga Food Allergen Lateral Flow Food reactivity data
Rapid Test ProⅡ for Shellfish

Rev. Mar. 2024

Lot number from #PSL002 -

Group	Ingredient	Test Result	Group	Ingredient	Test Result	Group	Ingredient	Test Result	
Cereals	Rice powder	-	Seeds	Peanut	-	Sea foods	Salmon	-	
	Polished rice	-		Roasted almond	-		Mackerel	-	
	Sprouted brown rice	-		Roasted cashew nut	-		Squid	+	
	Glutinous rice	-		Roasted macadamia nut	-		Octopus	+	
	Buckwheat	-		Roasted pistachio nut	-		Black tiger prawn	+	
	Japanese barnyard millet	-		Roasted walnut	-		Snow crab	+	
	Foxtail millet	-		Roasted pecan nut	-		Scallop adductor	+	
	Common millet	-		White sesame	-		Pollack roe	-	
	Wheat	-		Black sesame	-		Salmon roe	-	
	Rye	-	Green stuff	Corn	-	Meats	Beef(20%)	- *2	
	Barley	-		Tomato	-		Chicken	-	
	Malt	-		Onion	-		Pork(20%)	- *2	
	Oat	-		Spinach	-		Chicken egg	-	
	Adlay	- *1		Fruits	Orange	-	Eggs	Processed chicken egg	-
	Corn flour	-	Kiwi		-	Dried chicken egg		-	
	Amaranth	-	Apple		-	Dairy products		Milk	-
	Quinoa	-	Peach		-		Skim milk	-	
	Tubers・Starch	Sago palm powder	-		Fungus		Banana	-	Coffee・Cocoa・Tea
Yam		-	Coconut milk			-	Cocoa powder	-	
Potato starch		-	Mushroom	-		Green tea	+		
Beans	Soybean	-	Seaweeds	Roasted laver	-	Spices	Cumin	-	
	Red bean (azuki)	-		Hijiki (boiled)	-		Poppy seeds	-	
	Common bean	-		Brown seaweed (raw)	-		Black pepper	-	
	Concord pole	-		Kelp 200-fold extract	-		White pepper	-	
	Broad bean	-					Polysaccharide Thickener	Turmeric	-
								Low methoxyl pectin 50-fold extract	-
								Guar gum 1000-fold extract	-
								Gelatin	-
							Others		-

Test Result	
-	Negative
+	Positive

*1 The Test Solution may be difficult to flow on a test stick.
*2 The food extract was diluted five-fold (1:5 dilution) with a blank extraction buffer.

These samples were tested by following the instruction of Rapid Test ProⅡ (using Recommended Extraction Method), except for Kelp 200-fold extract,Polysaccharide and Thickener samples.

These results can be changed by the difference of kit's lot number or the sample's species or production area.