

Morinaga Food Allergen Lateral Flow  
Rapid Test Pro II for Buckwheat

Rev. Jan. 2022

Lot number from #PB001 -

Group	Ingredient	Test Result	Group	Ingredient	Test Result	Group	Ingredient	Test Result
Cereals	Rice powder	-	Sea foods	Salmon	-	Coffee・Cocoa・Tea	Coffee bean	-
	Polished rice	-		Mackerel	-		Cocoa powder	-
	Sprouted brown rice	-		Squid	-		Green tea	-
	Glutinous rice	-		Octopus	-	Seaweeds	Roasted laver	-
	Buckwheat			Black tiger prawn	-		Brown seaweed	-
	Japanese barnyard millet	-		Snow crab	-		Kelp 200-fold extract	+
	Foxtail millet	-		Scallop adductor	-	Spices	Cumin	-
	Common millet	+		Pollack roe	-		Poppy seeds	-
	Wheat	-		Salmon roe	-		Black pepper	-
	Rye	-	Meats	Beef	-		White pepper	-
	Barley	-		Chicken	-		Red pepper	-
	Malt	-		Pork	-		Wasabi (green horseradish)	-
	Oat	-	Eggs	Chicken egg	-		Ginger	-
	Adlay※	-		Processed chicken egg	-		Curry powder (spice mix)	-
	Corn	-		Dried chicken egg	-		Cinnamon	-
	Corn flour	-	Dairy products	Milk	-		Coriander	-
	Amaranth	-		Skim milk	-		Turmeric	-
	Quinoa	-	Green stuff	Tomato	-	Polysaccharide Thickener	Low methoxyl pectin 50-fold extract	-
Tubers・Starch	Sago palm powder	-		Onion	-		Guar gum 1000-fold extract	-
	Yam	-	Fungus	Spinach	-	Others	Gelatin	-
	Potato starch	-		Mushroom	-			
Beans	Soybean	-	Fruits	Orange	-			
	Red bean (azuki)	-		Kiwi	-			
	Common bean	-		Apple	-			
	Concord pole	-		Peach	-			
	Broad bean	-		Banana	-			
				Coconut milk	-			
Seeds	Peanut	-						
	Roasted almond	-						
	Roasted cashew nut	-						
	Roasted macadamia nut	-						
	Roasted pistachio nut	-						
	Roasted walnut	-						
	White sesame	-						
	Black sesame	-						

Test Result	
-	Negative
+	Positive (reaction positive)

N.T. : Not Tested  
※The Test Solution may be difficult to flow on a test stick.  
These samples were tested by following the instruction of Rapid Test Pro II for Buckwheat (using Recommended Extraction Method) except Polysaccharide and Thickener samples.

These results can be changed by the difference of kit's lot number or the sample's species or production area.

