

Morinaga Food Allergen Lateral Flow
Rapid Test Pro II for Crustacean

Rev. Feb. 2023

Lot number from #PR001 -

Group	Ingredient	Test Result	Group	Ingredient	Test Result	Group	Ingredient	Test Result
Cereals	Rice powder	-	Sea foods	Salmon	-	Coffee・Cocoa・Tea	Coffee bean	-
	Polished rice	-		Mackerel	-		Cocoa powder	-
	Sprouted brown rice	-		Squid	-		Green tea	-
	Glutinous rice	-		Octopus	-	Seaweeds	Roasted laver	-
	Buckwheat	-		Black tiger prawn	+ *1		Brown seaweed	-
	Japanese barnyard millet	-		Snow crab	+	Spices	Cumin	-
	Foxtail millet	-		Scallop adductor	-		Poppy seeds	-
	Common millet	-		Pollack roe	-		Black pepper	-
	Wheat	-		Salmon roe	-		White pepper	-
	Rye	-	Meats	Beef	-		Red pepper	-
	Barley	-		Chicken	-		Wasabi (green horseradish)	-
	Malt	-		Pork	-		Ginger	-
	Oat	-	Eggs	Chicken egg	-		Curry powder (spice mix)	-
	Adlay	- *2		Dried chicken egg	-		Cinnamon	-
	Corn	-	Dairy products	Milk	-		Coriander	-
	Corn flour	-		Skim milk	-		Turmeric	-
	Amaranth	-	Green stuff	Tomato	-	Polysaccharide Thickener	Low methoxyl pectin 50-fold extract	-
	Quinoa	-		Onion	-		Guar gum 1000-fold extract	-
Tubers・Starch	Sago palm powder	-		Spinach	-	Others	Gelatin	-
	Yam	-	Fungus	Mushroom	-			
	Potato starch	-		Orange	-			
Beans	Soybean	-	Fruits	Kiwi	-			
	Red bean (azuki)	-		Apple	-			
	Common bean	-		Peach	-			
	Concord pole	-		Banana	-			
	Broad bean	-		Coconut milk	-			
Seeds	Peanut	-						
	Roasted almond	-						
	Roasted cashew nut	-						
	Roasted macadamia nut	-						
	Roasted pistachio nut	-						
	Roasted walnut	-						
	White sesame	-						
	Black sesame	-						

Test Result	
-	Negative
+	Positive (reaction positive)

N.T. : Not Tested
*1 The test result may be negative if it is included in high concentrations.
*2 The Test Solution may be difficult to flow on a test stick.

These samples were tested by following the instruction of Rapid Test Pro II for Crustacean (using Recommended Extrction Method), except Kelp 200-fold extract, Polysaccharide and Thickener samples.

These results can be changed by the difference of kit's lot number or the sample's species or production area.