

Morinaga Food Allergen ELISA  
Soya ELISA Kit II

Rev. Nov. 2019

Lot number from #039( the last two or three digits )

| Group         | Ingredient               | Value<br>μg/g | Group          | Ingredient             | Value<br>μg/g | Group                    | Ingredient                            | Value<br>μg/g |
|---------------|--------------------------|---------------|----------------|------------------------|---------------|--------------------------|---------------------------------------|---------------|
| Cereals       | Rice powder              | <0.31         | Sea foods      | Salmon                 | <0.31         | Coffee・Cocoa・Tea         | Coffee bean                           | <0.31         |
|               | Polished rice            | <0.31         |                | Mackerel               | <0.31         |                          | Cocoa powder                          | <0.31         |
|               | Sprouted brown rice      | <0.31         |                | Tuna                   | <0.31         |                          | Green tea                             | <0.31         |
|               | Red rice                 | <0.31         |                | Dried bonito           | <0.31         |                          | Red tea                               | <0.31         |
|               | Black rice               | <0.31         |                | Squid                  | <0.31         | Seaweeds                 | Green laver                           | <0.31         |
|               | Glutinous rice           | <0.31         |                | Octopus                | <0.31         |                          | Roasted laver                         | <0.31         |
|               | Buckwheat                | <0.31         |                | Black tiger prawn      | <0.31         |                          | Brown seaweed                         | <0.31         |
|               | Japanese barnyard millet | <0.31         |                | Snow crab              | <0.31         |                          | Dried kelp                            | <0.31         |
|               | Foxtail millet           | <0.31         |                | Japanese spiny lobster | <0.31         | Spices                   | Rosemary                              | <0.31         |
|               | Common millet            | 0.35          |                | Lobster                | <0.31         |                          | Vanilla beans                         | <0.31         |
|               | Wheat                    | <0.31         | Meats          | Clam                   | <0.31         |                          | Cumin                                 | <0.31         |
|               | Rye                      | <0.31         |                | Abalone                | <0.31         |                          | Poppy seeds                           | <0.31         |
|               | Barley                   | <0.31         |                | Scallop adductor       | <0.31         |                          | Fennel                                | 0.64          |
|               | Malt                     | <0.31         |                | Pollack roe            | <0.31         |                          | Fapaniese peper                       | <0.31         |
|               | Oat                      | <0.31         |                | Salmon roe             | <0.31         |                          | Black pepper                          | <0.31         |
|               | Adlay                    | <0.31         |                | Raw beef               | <0.31         |                          | White pepper                          | <0.31         |
|               | Corn                     | <0.31         |                | Roasted beef           | <0.31         |                          | Cardamon                              | <0.31         |
|               | Corn flour               | <0.31         |                | Beef liver             | 2.19          |                          | Sage                                  | <0.31         |
|               | Amaranth                 | <0.31         |                | Raw chicken            | <0.31         |                          | Clove                                 | <0.31         |
|               | Quinoa                   | <0.31         |                | Roasted chicken        | <0.31         |                          | Red pepper                            | <0.31         |
| Tubers・Starch | Sago palm powder         | <0.31         |                | Chicken liver          | <0.31         |                          | Horseradish                           | <0.31         |
|               | Potato                   | <0.31         |                | Raw pork               | <0.31         |                          | Wasabi (green horseradish)            | <0.31         |
|               | Potato flake             | <0.31         |                | Roasted pork           | <0.31         |                          | Ginger                                | <0.31         |
|               | Yam                      | <0.31         |                | Pork liver             | <0.31         |                          | Curry powder (spice mix)              | 2.53          |
| Beans         | Soybean                  |               | Eggs           | Chicken egg            | <0.31         |                          | Nutmeg                                | <0.31         |
|               | Red bean (azuki)         | <0.31         |                | Processed chicken egg  | <0.31         |                          | Cinnamon                              | <0.31         |
|               | Common bean              | <0.31         |                | Dried chicken egg      | <0.31         |                          | Coriander                             | <0.31         |
|               | Green pea                | 0.59          | Dairy products | Milk                   | <0.31         |                          | Basil                                 | <0.31         |
|               | Runner bean              | <0.31         |                | Butter                 | <0.31         |                          | Thyme                                 | <0.31         |
|               | Kidney bean              | <0.31         |                | Skim milk              | <0.31         |                          | Laurel                                | <0.31         |
|               | Concord pole             | <0.31         | Green stuff    | Tomato                 | <0.31         |                          | Dried tangerine peel                  | <0.31         |
|               | Black bean               | >20           |                | Onion                  | <0.31         | Polysaccharide Thickener | Turmeric                              | <0.31         |
|               | Broad bean               | 1.07          |                | Spinach                | <0.31         |                          | Carrageenan (ι-type) 100-fold extract | <1.56         |
| Seeds         | Peanut                   | <0.31         | Fungus         | Garlic                 | <0.31         |                          | Carrageenan (κ-type) 50-fold extract  | <0.78         |
|               | Roasted almond           | <0.31         |                | Shiitake               | <0.31         |                          | Carrageenan (λ-type) 50-fold extract  | <0.78         |
|               | Roasted pecan nut        | <0.31         |                | Matsutake              | <0.31         |                          | Gellan gum 50-fold extract            | <0.78         |
|               | Roasted cashew nut       | <0.31         |                | Mushroom               | <0.31         |                          | Locust bean gum 200-fold extract      | 5.51          |
|               | Roasted macadamia nut    | <0.31         | Fruits         | Orange                 | <0.31         |                          | Low methoxyl pectin 50-fold extract   | <0.78         |
|               | Roasted pistachio nut    | <0.31         |                | Kiwi                   | <0.31         |                          | Xanthane gum 50-fold extract          | <0.78         |
|               | Roasted pistachio nut    | <0.31         |                | Apple                  | <0.31         |                          | Guar gum 1000-fold extract            | 128.4         |
|               | Roasted hazel nut        | <0.31         |                | Peach                  | <0.31         | Others                   | Gelatin                               | <0.31         |
|               | Roasted walnut           | 0.62          |                | Banana                 | <0.31         |                          | Yeast extract                         | <0.31         |
|               | White sesame             | <0.31         |                | Raisin                 | <0.31         |                          |                                       |               |
|               | Roasted white sesame     | <0.31         |                | Coconut powder         | <0.31         |                          |                                       |               |
|               | Black sesame             | <0.31         |                | Coconut milk           | <0.31         |                          |                                       |               |
|               | Roasted black sesame     | <0.31         |                |                        |               |                          |                                       |               |
|               | Ginkgo nut               | <0.31         |                |                        |               |                          |                                       |               |
|               | Pine nut                 | <0.31         |                |                        |               |                          |                                       |               |
|               | Matrimony vine seed      | 0.34          |                |                        |               |                          |                                       |               |
|               | Gardenia seed            | <0.31         |                |                        |               |                          |                                       |               |
|               | Hemp seed                | <0.31         |                |                        |               |                          |                                       |               |

N.T. : Not Tested

Reaction positive

※These samples were measured by following the instruction of this ELISA Kit.

※These values can be changed by the defference of kit's lot number or the sample's race or production area.