



FOOD ALLERGEN

Rapid Test Pro II



Features

- Lateral flow immunoassay kit
- Qualitative detection of food allergens
- Suitable for:
 - raw ingredients
 - highly processed foods
 - rinse water
 - environmental samples
- Ready to use-Extraction Solution
- Same sample extract can be analyzed in all Rapid Test Pro II kits
- All reagents and equipment (swab, tube, and disposable pipette) provided
- 24-month shelf life

Rapid Test Pro II	Specificity	Detection Limit for food	Detection Limit for swab (/100 cm ²)	Catalog No.
Egg	Ovalbumin	5 ppm Egg protein	1 µg Egg protein	M2261
Total Milk	Casein	4 ppm Casein	0.8 µg Casein	M2269
	β-lactoglobulin	0.5 ppm β-lactoglobulin	0.1 µg β-lactoglobulin	
Gluten	Gliadin	4 ppm Gluten	0.8 µg Gluten	M2263
Buckwheat	multiple buckwheat proteins	5 ppm Buckwheat protein	1 µg Buckwheat protein	M2264
Peanut	multiple peanut proteins	5 ppm Peanut protein	1 µg Peanut protein	M2265
Soya	β-conglycinin	5 ppm Soya protein	1 µg Soya protein	M2266
Crustacean	Tropomyosin	5 ppm Crustacean protein	1 µg Crustacean protein	M2267
Walnut	2S albumin	5 ppm Walnut protein	1 µg Walnut protein	M2268
Cashew	2S albumin	5 ppm Cashew protein	1 µg Cashew protein	M2260
Shellfish	Tropomyosin	5 ppm Shellfish protein	1 µg Shellfish protein	M2270


* ppm = µg protein / g of food

Procedure


Sample Preparation

Assay Procedure


Food Sample





Rinse water Sample



Swab Sample




Mix the Sample with Sample Extraction Solution

Dilute the Sample Extract 10-fold with Diluent





Apply 200 μ L (5 or 6 drops)



Vortex for mixing samples completely

For processed food

Expose in boiling water for 10 min and vortex





Result

Positive

Negative

{

}

For Egg and Shellfish

for Total Milk

Other kits

Kit components

Component	Amount	Component	Amount
Extraction Solution	10 packs	Polypropylene tube (L)	10
Diluent	1 bottle	Polypropylene tube (S)	10
Test stick	10 packs	Cotton swab	10 packs
Pipette (L)	10	Paper tube rack	1 piece
Pipette (S)	10	Insert	1 sheet



Morinaga BioScience, Inc.
 2-1-1 Shimosueyoshi,
 Tsurumi-ku, Yokohama-shi,
 230-8504
 Japan



<https://www.morinaga-biosci.com>



info_miobs_e@morinaga.co.jp



TEL +81-45-586-2515

FAX +81-45-586-2517