

Morinaga Institute of Biological Science, Inc.

SACHIURA 2-1-16, KANAZAWA-KU, YOKOHAMA, 2360003, JAPAN TEL:+81-45-791-7674 E-mail: info_miobs_e@morinaga.co.jp

TEL:+81-45-791-7674 E-mail: info_miobs_e@morinaga.co.jp https://www.miobs-e.com

Announcement of product revision

October. 2022

Dear Customers,

Morinaga Institute of Biological Science, Inc. now announces that a series of Rapid Test Pro II has been revised, as shown in the tables below. However, there is no change in the Limit of detection and the procedures. The composition of the Diluent contained in the Rapid Test Pro II for Egg, Total Milk, Gluten, and Peanut differs from the ones of the former products. So please do not mix-up the former Diluent.

Before revision

Product	Cat. No.	Lot #	Diluent	Positive Line
Rapid Test Pro II for Egg	M2261	∼ PE022	In 3 products, Diluent is the same composition (Red cap bottle)	Red
Rapid Test Pro II for Casein	M2262	~ PC022		
Rapid Test Pro II for Gluten	M2263	∼ PG022		
Rapid Test Pro II for Peanut	M2265	∼ PP020	In 2 products, Diluent is the same composition (Blue cap bottle)	
Rapid Test Pro II for Buckwheat	M2264	∼ PB021		

After revision

Product	Cat. No.	Lot #	Diluent	Positive Line
Rapid Test Pro II for Egg	M2261	PE023 \sim	In 4 products, Diluent is the same composition (Clear cap bottle)	Black
Rapid Test Pro II for Total Milk	M2269	PM001 ∼		Red & Black
Rapid Test Pro II for Gluten	M2263	PG023 ∼		Red
Rapid Test Pro II for Peanut	M2265	PP021 ∼		
Rapid Test Pro II for Buckwheat	M2264	PB022 ∼	Only for Buckwheat kit (Blue cap bottle)	

^{*}Since Rapid Test Pro II for Total Milk, M2269, has been launched, Rapid Test Pro II for Casein, M2262, is discontinued.

If you have any questions, please contact info miobs e@morinaga.co.jp.

Sincerely yours,

Sales team, Morinaga Institute of Biological Science, Inc.