

Morinaga Food Allergen ELISA Kit II Food reactivity data
Peanut ELISA Kit II

Rev. Sep-25

Lot number from			FPE090	~	now			
Group	Ingredient	Value μg/g	Group	Ingredient	Value μg/g	Group	Ingredient	Value μg/g
Cereals	Rice powder	<0.31	Green stuff	Corn	<0.31	Dairy products	Milk	<0.31
	Polished rice	<0.31		Tomato	<0.31		Butter	<0.31
	Sprouted brown rice	<0.31		Onion	<0.31		Skim milk	<0.31
	Red rice	<0.31		Spinach	<0.31	Coffee・Cocoa・Tea	Coffee bean	<0.31
	Black rice	<0.31		Garlic	<0.31		Cocoa powder	<0.31
	Glutinous rice	<0.31	Fruits	Orange	<0.31		Green tea	<0.31
	Buckwheat	<0.31		Kiwi	<0.31	Spices	Black tea	<0.31
	Japanese barnyard millet	0.35		Apple	<0.31		Rosemary	<0.31
	Foxtail millet	0.54		Peach	<0.31		Vanilla beans	<0.31
	Common millet	0.75		Banana	<0.31		Cumin	<0.31
	Wheat	<0.31		Raisin	<0.31		Poppy seeds	<0.31
	Rye	<0.31		Coconut powder	<0.31		Fennel	<0.31
	Barley	<0.31		Coconut milk	<0.31		Japanese pepper	<0.31
	Malt	<0.31	Fungus	Shiitake	<0.31		Black pepper	<0.31
	Oat	<0.31		Matsutake	<0.31		White pepper	<0.31
	Adlay	<0.31		Mushroom	<0.31		Cardamon	<0.31
	Corn flour	<0.31	Seaweeds	Green laver	<0.31		Sage	<0.31
	Amaranth	<0.31		Roasted laver	<0.31		Clove	<0.31
	Quinoa	<0.31		Brown seaweed	<0.31		Red pepper	<0.31
	Sago palm powder	<0.31		Dried kelp	<0.31		Horseradish	<0.31
Tubers・Starch	Potato	<0.31	Sea foods	Salmon	<0.31		Wasabi (green horseradish)	<0.31
	Potato flake	<0.31		Mackerel	<0.31		Ginger	<0.31
	Yam	<0.31		Tuna	<0.31		Curry powder (spice mix)	<0.31
				Dried bonito	<0.31		Nutmeg	<0.31
Beans	Soybean	<0.31		Squid	<0.31	Polysaccharide Thickener	Cinnamon	<0.31
	Red bean (azuki)	<0.31		Octopus	<0.31		Coriander	<0.31
	Common bean	<0.31		Black tiger prawn	<0.31		Basil	<0.31
	Green pea	<0.31		Snow crab	<0.31		Thyme	<0.31
	Runner bean	<0.31		Japanese spiny lobster	<0.31		Laurel	<0.31
	Kidney bean	<0.31		Lobster	<0.31		Dried tangerine peel	<0.31
	Concord pole	<0.31		Clam	<0.31		Turmeric	<0.31
	Black bean	<0.31		Abalone	<0.31		Carrageenan (ι-type) 100-fold extract	<1.56
	Broad bean	<0.31		Scallop adductor	<0.31		Carrageenan (κ-type) 50-fold extract	<0.78
Seeds	Peanut	>20		Pollack roe	<0.31		Carrageenan (λ-type) 50-fold extract	<0.78
	Roasted almond	<0.31	Meats	Salmon roe	<0.31		Gellan gum 50-fold extract	<0.78
	Roasted pecan nut	<0.31		Raw beef	<0.31		Locust bean gum 200-fold extract	<3.12
	Roasted cashew nut	<0.31		Roasted beef	<0.31		Low methoxyl pectin 50-fold extract	<0.78
	Roasted macadamia nut	1.2		Beef liver	<0.31		Xanthane gum 50-fold extract	<0.78
	Roasted pistachio nut	<0.31		Raw chicken	<0.31	Others	Guar gum 1000-fold extract	<15.6
	Roasted hazel nut	<0.31		Roasted chicken	<0.31		Gelatin	<0.31
	Roasted walnut	<0.31		Chicken liver	<0.31		Yeast extract	<0.31
	White sesame	<0.31		Raw pork	<0.31			
	Roasted white sesame	<0.31		Roasted pork	<0.31			
	Black sesame	<0.31		Pork liver	<0.31			
	Roasted black sesame	<0.31	Eggs	Chicken egg	<0.31			
	Ginkgo nut	<0.31		Processed chicken egg	<0.31			
	Pine nut	<0.31		Dried chicken egg	<0.31			
	Matrimony vine seed	<0.31						
	Gardenia seed	<0.31						
	Hemp seed	<0.31						

Reaction positive

※These samples were measured by following the instruction of this ELISA Kit.

※These values can be changed by the difference of kit's lot number or the sample's race or production area.