



Morinaga BioScience, Inc.
2-1-1 SHIMOSUEYOSHI, TSURUMI-KU,
YOKOHAMA, 230 8504, JAPAN
TEL: +81-45-586-2515
E-mail: info_miobs_e@morinaga.co.jp

Announcement of the 1 ppm LOD Testing Protocol in Rapid Test Pro II

October, 2026

To Our Valued Customers,

We are writing to announce two important updates regarding the Rapid Test Pro II product line: the addition of the 30-minute method for select products, and the transition from paper instruction manuals to digital manuals across the entire series.

1. Addition of the 30-minute method (1 ppm Limit of Detection Testing)

For three products in the Rapid Test Pro II series (Egg, Total Milk, Gluten), the 30-minute method has been added to enable detection at 1 ppm.

The sample preparation and handling procedure remains identical to the current method. By simply adjusting the incubation time after applying the Test Solution, users can now select either 5 ppm (15 min incubation) or 1 ppm (30 min incubation) detection sensitivity according to their testing requirements. Please do not interpret the test results between 15 and 30 min. (e.g. at 20 min or 25 min).

2. Digitalization of Kit Instructions

As part of our eco-friendly initiatives, we are transitioning from paper manuals to digital instructions for all nine products in the Rapid Test Pro II series. A paper insert containing a 2D code will be provided inside each kit box. Please scan the 2D code below with your mobile device or visit our official website to view or download the latest kit instruction.



These two updates will apply starting with the lot numbers listed in the table below.



Morinaga BioScience, Inc.
 2-1-1 SHIMOSUEYOSHI, TSURUMI-KU,
 YOKOHAMA, 230 8504, JAPAN
 TEL: +81-45-586-2515
 E-mail: info_miobs_e@morinaga.co.jp

Product Name	Cat. No.	30-Minute Method (1 ppm testing)	Digital Manual Transition	Effective Lot No.
Rapid Test Pro II for Egg	M2261	Available (5 ppm testing / 1 ppm testing)	Applied (2D code)	PEG042
Rapid Test Pro II for Gluten	M2263			PGT040
Rapid Test Pro II for Total Milk	M2269			PTM022
Rapid Test Pro II for Cashew	M2260	PCW005		
Rapid Test Pro II for Buckwheat	M2264	PBW038		
Rapid Test Pro II for Peanut	M2265	PPE035		
Rapid Test Pro II for Soya	M2266	PSY028		
Rapid Test Pro II for Crustacean	M2267	PCR014		
Rapid Test Pro II for Walnut	M2268	PWL017		
Rapid Test Pro II for Shellfish	M2270			PSL008

If you have any questions, please contact us at info_miobs_e@morinaga.co.jp.

Faithfully yours,
 Sales team
 Morinaga BioScience, Inc.